EFM: Journey Through Time

Emily Hamilton MDCM July 17, 2015

Contents

Introduction	3
Evolution of Electronic Fetal Monitoring	4
Sensor Phase	4
Physiology Phase	5
Guideline Phase	6
Performance of EFM	6
Future Paths	10
Conclusions	11
References	11
Appendix of EFM Classification Methods	14
1. EFM scoring system based on Krebs et al 1979	15
2. Classification based on FIGO methods in 1987	16
3. Classification and management based on guidelines by ACOG 2010	17
4. Tracing classification and management based on NICE guidelines	18
5. Classification and management by Parer et al 2006	19
6. Management algorithm for Category II tracings by Clark et al 2012.	20

Introduction

Controversy is no stranger to electronic fetal monitoring (EFM). Few technologies in medicine can claim association with such a wide range of professional reaction. At different times EFM technology has been inflated, berated, debated and sometimes obfuscated. Despite controversy, it remains a mainstay of intrapartum care, suggesting that clinicians find its benefits outweigh its disadvantages.

Clinical benefit is supported by many encouraging reports show falling rates of intrapartum-related neonatal encephalopathy (NE). A myriad of factors influence complicated outcomes such as NE and it would be incorrect to attribute the improvement to any single one. That said, several large jurisdictions with population-wide statistics have shown steady declines, often reaching reductions in the range of 40% per decade. Figure 1 shows declining numbers of births with intrapartum-related neonatal encephalopathy by decade in several regions. (7)




Figure 1. Regional trends in the numbers of births with intrapartum-related neonatal encephalopathy

Improvement in outcome is not likely to be related simply to better recording of the fetal heart rate. EFM monitors have changed little in the past thirty years. In contrast, clinical practices today are very different from those in the mid 1980's when the largest randomized clinical trial (RCT) on EFM was reported. During that study, perinatal hypoxic death or newborn seizures occurred at an astounding rate of approximately 1 in every 225 births and primary cesarean rates were around 2.4%. (8) Today EFM is used in the vast majority of births in hospitals. In high-

income countries, the rate of NE is around 1.5 per 1000 births. Preventable intrapartum stillbirths are almost eliminated and cesarean rates often exceed 30%.

Over the last few decades we have come to realize that that human factors and system failures play a substantial role in adverse outcomes across all of branches medicine. Human actions, such as delayed recognition of tracing abnormality or delayed intervention, are reported to have occurred in approximately half of birth-related asphyxia injuries. Drawing from aviation and military experience, we have adopted their models to build less error-prone health care systems. Consequently policies and procedures to redress the reasons underlying human error are now widespread. In middle- and low-income countries, improving the national socioeconomic status and access to basic and safe healthcare is strongly associated with improved health outcomes.

In order to better understand the relationship between the EFM and the improvement in NE it is helpful to review the timeline of advances in fetal monitoring.

Evolution of Electronic Fetal Monitoring

The invention of the stethoscope is generally attributed to René Laennec in 1816. He listened to adult heart tones using a roll of paper and later fashioned a wooden tube-like instrument. "Ear trumpets" (flared funnel shaped devices used to concentrate sound waves) had been used as hearing aids for more than a century. His internist friend, Jean Alexandre Lejumeau, Vicomte de Kergaradec, tried this device to listen to the "noises" of amniotic fluid but instead he heard fetal heart tones. His questioning merits translation: "From the changes occurring in the strength and rate of fetal heart beats, wouldn't it be possible to know about the status of health or sickness of the fetus?" (19, 20) How prescient considering that we ask the same question nearly 200 years later!

Despite the challenges of disseminating medical information in that era, the significance of fetal heart tones and heart rate as an indicator of fetal status was recognized relatively rapidly. By 1833 Dublin physician Evory Kennedy published a monograph on auscultation of the fetal heart, commenting on the ominous heart rate pattern of "slowness of its return when a contraction is passing on," the effects of head or cord compression on heart rate, and the significance of meconium stained amniotic fluid. (20, 21) The first stethoscope designed specifically for fetal heart tones was created 79 years later by Adolphe Pinard in 1895. Another half century passed before reports appeared describing electronic methods to detect the fetal heart rate. Figure 2 shows a timeline of three phases and key events in the evolution in electronic fetal monitoring.

Sensor Phase

Spanning two decades, this first phase saw the creation of the basic electronic sensors and equipment to measure the fetal heart rate (FHR) and contractions. By 1968 the first commercial fetal monitor device was released using phonocardiography and external tocography.

Phonocardiography quickly gave way to direct scalp electrode for the measurement of FHR followed by the appearance of ultrasound-based sensors (22)




Figure 2. Historical phases in the evolution of electronic fetal monitoring

Physiology Phase

During the next three decades, our understanding of human fetal physiology grew extensively. With measurement now possible it became feasible to observe fetal heart rate patterns during pregnancy and in labor. Using animal experimentation with pharmacological blocking agents, vessel cannulation/occlusion and nerve transection studies researchers elucidated how the fetal heart was regulated. Retrospective analysis of tracing in labors with adverse outcomes or other "natural experiments" like anencephaly or severe fetal anemia revealed other insights. (Important cardiac control mechanisms are discussed in more detail in the document, *EFM Basics: Physiology.*) During the 30-year physiology phase, over 8,000 publications are found in PubMed using the search terms fetal heart rate. An outstanding summary and organization of EFM literature was compiled by the Royal College of Obstetricians and Gynaecologists in the UK and remains an excellent resource today (23)

Guideline Phase

As information became available from diverse research teams and electronic fetal heart rate monitoring adoption rose, an assortment of new terms followed. By 1971 the first international conference on monitoring of the fetal heat convened and developed consensus on related terminology. Unfortunately, their conclusions were never published. Over the years terminology was revisited by several professional groups. (24-28) Measuring the FHR, labeling patterns, and understanding the relevant physiology did not in themselves improve outcomes. However, they were the foundation upon which all guidelines directing clinical management were built.

Clinical guidelines for EFM management generally have two components:

- A graded classification defines the degree of tracing abnormality which in turn ordains the nature and urgency of intervention.
- A variety of EFM scoring systems were developed for both the antepartum and intrapartum periods. The newly introduced term *reactive* antepartum tracing gave rise to no less than 21 different definitions of a reactive test!⁽²⁹⁾ Early methods for the more difficult task of intrapartum classification used scoring systems. Points were awarded for various features with a maximum possible score of 10. Higher scores were more reassuring than lower scores. In practice scores were generally grouped into three levels so in effect 3-level classification systems have been is use for about 40 years.

Early classification systems published by Krebs et al in 1979 and FIGO in 1987 are summarized Appendix 1. $^{(30, 31)}$

Four modern classification systems place more reliance upon minimal baseline variability as a deciding factor, specify more details on the size, number and type of decelerations that are of concern and recommend clinical actions. For ease of comparison, they are also summarized in Appendix 1. (32-35)

Performance of EFM

The clinical goals of EFM are to identify fetuses with increased risk of hypoxic injury so that intervention can avoid adverse outcome without also causing excessive numbers of unnecessary interventions. Thus, it is important to measure how often the classification method detected as abnormal the tracings of babies who showed hypoxic injury and how often it did the same in babies who were born without a problem. These measures are the hallmarks used to assess the performance of a diagnostic method. However the goal of EFM is *prevention* of illness and not diagnosis. It is not straightforward to measure the performance of a prevention technique such as EFM for two important reasons.

The first reason is related to "intervention paradox". The fundamental measures of performance of a diagnostic are sensitivity (% of unhealthy patients with a positive test) and false positive rates (% of healthy patients with a positive test). With EFM we have an intervention paradox. When the EFM is positive, intervention can prevent the illness from

occurring. The intervention paradox is that a *positive* test is now accompanied by a *healthy* outcome due to successful clinical intervention. We rarely know with certainty in this situation if a bad outcome was truly averted or if the EFM-based indication was just a false positive. Unless there is a test that can indicate that a bad outcome was impending, intervention paradox will cause us to underestimate sensitivity and overestimate false positive rates. Thus we should keep in mind that sensitivity and specificity are very conservative estimates of the value of FFM.


Intervention paradox should not prevent us from measuring sensitivity and specificity. These fundamental performance measures remain useful especially when comparing EFM classification systems done on the same dataset because they are subjected to the same limitations. Both high sensitivity and high specificity are desirable for different reasons. High sensitivity is desirable because hypoxic injury can have devastating long-term consequences. High specificity is desirable because normal outcome is far more common than hypoxic injury and many interventions based on a high false positive rate will create a major health care burden with no benefit. Summing sensitivity and specificity is a simple way of combining both measures with equal weight in order to compare the classification methods. A perfect test would have a score of 200; a test that was no better than chance would have a score of 100.


Table 1 shows a summary of performance for a variety of classification methods. (36-40) The sensitivity and false positive rates are shown graphically in Figure 3. The studies are not comparable in terms of adverse outcome studied. Some examined HIE, others used various levels of acidemia. Some used visual analysis and others used automated methods. Nevertheless the results show an interesting pattern.

Table 1. Summary of reports describing classification performance

Source	Outcome Number	Control Number	Classification Method	Sensitivity	Specificity	False positive	Sum of Sensitivity and Specificity
Spencer et al 1997	HIE all grades	No HIE	FIGO levels suspicious or abnormal	89%	48%	52%	137%
	38	35					
Low et al 1999	UA base deficit > 16 mmol/L	UA base deficit <8mmol/L					
	71	71	Absent baseline variability > 10 min	17%	98%	2%	115%
			Minimal variability > 20 min AND late or	46%	89%	11%	135%
			prolonged decelerations Minimal variability >20 min OR late and/or prolonged decelerations	75%	57%	43%	132%
Elliott et al 2006	HIE all grades	UA base deficit <8 mmol/L	Parer and Ikeda 5 levels				
	60	2132	Red	8%	98.3%	1.7%	107%
			Orange and above	52%	81.9%	18.1%	134%
			Yellow and above	88%	40.1%	59.9%	128%
			Blue and above	98%	13.1%	86.9%	111%
Elliott et al 2006	UA base deficit > 12 mmol/L	UA base deficit <8 mmol/L	Parer and Ikeda 5 levels				
	280	2132	Red	3%	98.3%	1.7%	101%
			Orange and above	30%	81.9%	18.1%	112%
			Yellow and above	73%	40.1%	59.9%	113%
			Blue and above	94%	13.1%	86.9%	107%
Cahill et al 2012	UA pH<=7.10	UA pH>7.10	ACOG 3 Categories				
	57	5331	Category II	100%	2.4%	97.6%	102%
			Category III	0%	99.9%	0.1%	100%
Clark et al 2015	UA base deficit > 12 mmol/L	UA base deficit <8 mmol/L	9branches based on EFM and labor status				
	120	120	Any Intervention Recommendation	45.8%	81.7%	18.3%	128%

Figure 3 demonstrates the limitation of rule-based classification methods using EFM and EFM classical features. High sensitivity can be achieved only at the expense of a high false positive rate. The general relationship of sensitivity and specificity is shown by the solid line. Any new technique with a true advance should show performance levels much higher and to the left of this line. The best performing method is represented by the spot that is closest to the upper left corner of the graph.





The second reason measuring the performance of EFM is challenging relates to the use of clinical outcome as a measurement. Clinical outcome results from the cumulative effect of several steps namely beginning with signal acquisition, followed by diagnosis, and ending with timely and effective intervention. A deficiency at any step can negate all the benefit of a previous step. Likewise performance at an early step can hamper all subsequent steps. Several studies have reported that human actions, such as delayed recognition of tracing abnormality and/or delayed intervention, were present in large percentages of cases with asphyxial injuries. (10-14) Thus, clinical outcome measures reflect the entire process instead of any single step. Reflect upon the similarity of the basic electronic fetal monitor and the classifications systems used in the 1980s to their counterparts in place today. Contrast this with the rate of perinatal death or HIE of 1 in 225 in the 1980s to 1.5 per 1000 commonly seen today. Something happened to improve outcomes.

The patient safety movement has helped us understand the causes of error-prone health care systems. Some mitigation is directed at system vulnerabilities, such as legislation to limit working hours, recommendations on staff-to-patient ratios, setting standards for availability of obstetricians and operating room facilities, simulation training for emergency procedures and formal feedback on performance to clinicians. Other actions are very specific to electronic fetal monitoring: Standardized nomenclature; graded classifications of abnormality; and formal guidelines for clinical management. Finally, improved understanding about the clinical significance of associating some fetal heart rate patterns with outcome helps clinicians respond

better. For example defining the correlation between rates of fetal death or HIE and the interval between persistent bradycardia and delivery is crucial to establishing desirable response times. (41)

Another advance has been the growing realization that small changes in critical areas can have a major impact. Simple checklists, limited to a few top items have been associated with dramatic improvements across diverse medical specialties. Some of the most notable checklist achievements are highlighted in Figure 4.⁽⁴²⁾




Figure 4. Falling complication rates associated with clinical checklist use.

In obstetrics, compliance with a simple checklist of six EFM related conditions to trigger discontinuation of oxytocin was associated with significant clinical benefit. In a study spanning four years and 14,398 term inductions, outcomes were compared in patients with and without checklist compliance. (43) Compliance was associated with fewer NICU admissions (4.4% vs 2.9%) and a lower primary cesarean rate (18.8% vs 15.8%). Given all of the limitations of EFM, usage consistent with management guidelines was associated with improved outcomes.

Future Paths

Today clinicians rely largely upon visual inspection of electronic fetal monitoring (EFM) tracings to assess fetal tolerance to labor. Rule-based classification methods are limited since they require a tradeoff between sensitivity and false positive rates. The EFM features in the rules are constrained to what the human eye can recognize and measure. Human assessment is inconsistent, especially when assessments are carried over long periods of time in the presence of clinician fatigue or distractions. Low incidence-high consequence medical problems are among the most challenging for clinicians, especially when false positives are common. A FHR aberration is almost always associated with a good outcome from the perspective of a front-line clinician. Considering all of these impediments it is quite remarkable that HIE levels have fallen, albeit at the cost of high cesarean rates.

Active exploration in three areas is likely to advance the quest to reduce birth-related injury.

- Organizational and human behavior science will help us design and sustain higher
 quality healthcare systems. However, even a perfectly functioning labor & delivery unit
 cannot overcome the basic limitation of current EFM guidelines to identify and
 intervene on behalf of fetuses with increased risk of hypoxic injury without excessive
 numbers of unnecessary interventions.
- New sensors could provide a better or more direct measure of fetal cerebral state with respect to impending hypoxic injury. Experience with fetal oxygen saturation techniques and ST segment analysis of the fetal ECG have been disappointing to date.
- 3. Better analysis of existing EFM signals using approaches that can measure components that are not readily measurable or visible to the human eye.

Conclusions

We are entering a new phase in the evolution of fetal surveillance. Historically EFM research was stymied due to our inability to access and analyze the large data sets. We have entered the "big data" era in obstetrics. Hospitals store vast amounts of clinical data and digital EFM tracings. We can now analyze these digital EFM signals directly, measure standard EFM features as well as components and relationships that are not readily visible to humans. As these impediments are overcome we are seeing a resurgence of research producing a better understanding of which EFM characteristics are truly predictive of neonatal depression or metabolic acidosis and a search for new sensors. (44-46)

We all look forward to the next phase in the evolution of fetal surveillance where we will have better healthcare delivery practices and better technology to help with the difficult challenge of assessing the fetal status during labor.

References

- 1. Cyr RM, Usher RH, McLean FH. Changing patterns of birth asphyxia and trauma over 20 years. Am J Obstet Gynecol. 1984 Mar 1;148(5):490-8.
- 2. Hull J, Dodd KL. Falling incidence of hypoxic-ischaemic encephalopathy in term infants. Br J Obstet Gynaecol. 1992 May;99(5):386-91.
- 3. Smith J, Wells L, Dodd K. The continuing fall in incidence of hypoxic-ischaemic encephalopathy in term infants. BJOG. 2000 Apr;107(4):461-6.
- 4. Bell R, Glinianaia SV, Rankin J, Wright C, Pearce MS, Parker L. Changing patterns of perinatal death, 1982-2000: a retrospective cohort study. Arch Dis Child Fetal Neonatal Ed. 2004 Nov;89(6):F531-6.
- 5. Becher JC, Stenson BJ, Lyon AJ. Is intrapartum asphyxia preventable? BJOG. 2007 Nov;114(11):1442-4. Epub 2007 Sep 17.

- 6. Pasupathy D, Wood AM, Pell JP, Fleming M, Smith GC. Rates of and factors associated with delivery-related perinatal death among term infants in Scotland. JAMA. 2009 Aug 12;302(6):660-8.
- Lee ACC, Kozulki N, Blencowe H et al. Intrapartum-related neonatal encephalopathy incidence and impairment at regional and global levels for 2010 with trends from 1990. PedResearch 2013;74: 50-72.
- 8. MacDonald D, Grant A, Sheridan-Pereira M, Boylan P, Chalmers I. The Dublin randomized controlled trial of intrapartum fetal heart rate monitoring. Am J Obstet Gynecol. 1985 Jul 1;152(5):524-39.
- 9. Reason JT. (2001) Understanding Adverse Events: the Hyman Factor in C. Vincent (ed) Clinical Risk Management. London, UK BMJ Books
- 10. Clark SL, Belfort MA, Dildy GA, Meyers JA. Reducing obstetric litigation through alterations in practice patterns. Obstet Gynecol. 2008 Dec;112(6):1279-83.
- 11. Joint Commission on Accreditation of Healthcare Organizations, USA. Preventing infant death and injury during delivery. Sentinel Event Alert. 2004;(30):1-3.
- 12. Ransom SB, Studdert DM, Dombrowski MP, Mello MM, Brennan TA. Reduced medicolegal risk by compliance with obstetric clinical pathways: a case—control study. Obstet Gynecol. 2003;101(4):751-5.
- 13. 1Draper ES, Kurinczuk JJ, Lamming CR, Clarke M, James D, Field D. A confidential enquiry into cases of neonatal encephalopathy. Arch Dis Child Fetal Neonatal Ed. 2002;87(3):F176-80.
- 14. Saphier CJ, Thomas EJ, Studdert D, Brennan TA, Acker D. Applying no-fault compensation criteria to obstetric malpractice claims. Prim Care Update Ob Gyns. 1998;5(4):208-9.
- 15. Pettker CM, Thung SF, Raab CA, Donohue KP, Copel JA, Lockwood CJ, Funai EF. A comprehensive obstetrics patient safety program improves safety climate and culture. Am J Obstet Gynecol. 2011 Mar;204(3):216.e1-6.
- 16. Clark SL, Meyers JA, Frye DK, Perlin JA. Patient safety in obstetrics--the Hospital Corporation of America experience. Am J Obstet Gynecol. 2011 Apr;204(4):283-7.
- 17. Bohmer RM. The four habits of high-value health care organizations. N Engl J Med. 2011 Dec 1;365(22):2045-7.
- 18. Knox GE, Simpson KR. Perinatal high reliability. Am J Obstet Gynecol. 2011 May;204(5):373-7.
- 19. Kergaradec JA Lejumeau de (1822) Mémoire sur l'ausculation appliquée à l'etude de la grossesse. Communication Académie Médecine Paris 26/12/1821 p 1-37 Paris:Mequignon, Marves.
- 20. Sureau C. Historical persepecitives: forgotton past, unpredictable future. Ballière's Clinical Obstetrics and Gynecology 1996;10:2:167-84.
- 21. Kennedy E. Observations of Obstetrical Auscultation. Dublin: Hodges and Smith
- 22. Freeman RK , Garite TJ, Nageotte MP, Miller LA. Fetal Heart Rate Monitoring. 4th ed. Philadelphia, PA: Lippincott, Williams and Wilkins;2012
- 23. Royal College of Obstetricians and Gynecologists. Electronic fetal monitoring: the use and interpretation of cardiotocography in intrapartum fetal surveillance. Evidence-based guideline no. 8. Available at: http://guidance.nice.org.uk/CGC. Accessed June 30, 2006.
- 24. Predictors of intrapartum fetal distress: the role of electronic fetal monitoring. Report of the National Institute of Child Health and Human Development Consensus Development Task Force. Zuspan FP, Quilligan EJ, Iams JD, van Geijn HP. Am J Obstet Gynecol. 1979 Oct 1;135(3):287-91.
- 25. Electronic fetal heart rate monitoring: research guidelines for interpretation. National Institute of Child Health and Human Development Research Planning Workshop. Am J Obstet Gynecol. 1997 Dec;177(6):1385-90.
- 26. International Federation of Gynecology and Obstetrics. Guidelines for the use of fetal monitoring. Int J Gynaecol Obstet 1987;25:159-67.

- 27. Society of Obstetricians and Gynecologists of Canada. Fetal Health Surveillance in Labour. Ottawa: SOGC Policy Statement no. 41.
- 28. Fetal health surveillance in labour. Liston R, Crane J, Hamilton E, Hughes O, Kuling S, MacKinnon C, McNamara H, Milne K, Richardson B, Trépanie MJ; Working Group on Fetal Health Surveillance in Labor, Executive and Council, Maternal-Fetal Medicine Committee, Clinical Practice Guidline Committee, and ALARM Committee, Society of Obstetricians and Gynaecologists Canada; Canadian Medical Protection Association. J Obstet Gynaecol Can. 2002 Mar;24(3):250-76.
- 29. Devoe LD, Castillo RA, Sherline DM. The nonstress test as a diagnostic test: a critical reappraisal. Am J Obstet Gynecol. 1985 Aug 15;152(8):1047-53.
- 30. Krebs HB, Petres RE, Dunn LJ, Jordaan HV, Segreti A. Intrapartum fetal heart rate monitoring. I. Classification and prognosis of fetal heart rate patterns. Am J Obstet Gynecol. 1979 Apr 1;133(7):762-72.
- 31. FIGO News. Guidelines for the use of fetal monitoring, Int J Gynaecol Obstet 1987; 25: 159–167.
- 32. American College of Obstetricians and Gynecologists. Practice bulletin no. 116: Management of intrapartum fetal heart rate tracings. Obstet Gynecol. 2010 Nov;116(5):1232-40.
- 33. Nice Guidelines for fetal monitoring Available at:
 http://www.nice.org.uk/guidance/cg190/chapter/1-recommendations#monitoring-during-labour
 accessed July 10, 2015
- 34. Parer JT, Ikeda T. A framework for standardized management of intrapartum fetal heart rate patterns. Am J Obstet Gynecol. 2007;197:26.e1-6.
- 35. Clark SL, Nageotte MP, Garite TJ, Freeman RK, Miller DA, Simpson KR, Belfort MA, Dildy GA, Parer JT, Berkowitz RL, D'Alton M, Rouse DJ, Gilstrap LC, Vintzileos AM, van Dorsten JP, Boehm FH, Miller LA, Hankins GD. Intrapartum management of category II fetal heart rate tracings: towards standardization of care. Am J Obstet Gynecol. 2013 Aug;209(2):89-97.
- 36. Spencer JA, Badawi N, Burton P, Keogh J, Pemberton P, Stanley F. The intrapartum CTG prior to neonatal encephalopathy at term: a case-control study. Br J Obstet Gynaecol. 1997 Jan;104(1):25-8.
- 37. Low JA, Victory R, Derrick EJ. Predictive value of electronic fetal monitoring for intrapartum fetal asphyxia with metabolic acidosis. Obstet Gynecol. 1999 Feb;93(2):285-91.
- 38. Elliott C, Warrick PA, Graham E, and Hamilton EF. Graded classification of fetal heart rate tracings: Association with neonatal metabolic acidosis and neurologic morbidity. Am J Obstet Gynecol. 2010 Mar;202(3):258.e1-8.
- 39. Cahill AG, Roehl KA, Odibo AO, Macones GA. Association and prediction of neonatal acidemia. Am J Obstet Gynecol. 2012 Sep;207(3):206. e1-8.
- 40. Clark et al under review
- 41. Leung TY, Chung PW, Rogers MS, Sahota DS, Lao TT, Hung Chung TK. Urgent cesarean delivery for fetal bradycardia. Obstet Gynecol. 2009 Nov;114(5):1023-8.
- 42. Gawande, A. A lifesaving checklist. The New Yorker, December 10, 2007.
- 43. Clark SL, Meyers JA, Frye DK, Garthwaite T, Lee AJ, Perlin JB. Recognition and response to electronic fetal heart rate patterns: impact on newborn outcomes and primary cesarean delivery rate in women undergoing induction of labor. Am J Obstet Gynecol. 2015 Apr;212(4):494.e1-6.
- 44. Warrick PA, Hamilton EF. Discrimination of normal and at-risk populations from fetal heart rate variability. Computing in Cardiology 2014;41:1001-4.
- 45. Chudacek V, Anden J, Mallat S et al. Scattering transform for intrapartum fetal heart rate variability fractal analysis: a case-control study. IEEE Transactions on Biomedical Engineering 2014;61:1100-8.
- 46. Warrick PA, Hamilton EF, Precup D, Kearney RE. Classification of normal and hypoxic fetuses from systems modeling of intrapartum cardiotocography. IEEE Trans Biomed Eng. 2010 Apr;57(4):771-9.

Appendix of EFM Classification Methods


1. EFM scoring system based on Krebs et al 1979

Krebs Score							
FHR Feature	2	1	0				
Baseline FHR (bpm)	120-160	100-119	< 100				
, , ,		161-180	>180				
Variability							
Amplitude (bpm)	6-25	3-5 or > 25	<3				
Frequency	>6	3-6	<3				
oscillations/min							
Accelerations	>4	1-4	0				
Decelerations	None	Moderate Variable	Late				
	Early		Severe Variable Atypical variable				

2. Classification based on FIGO methods in 1987

	FIGO Classification					
FHR Feature	Normal	Suspicious	Pathological			
Baseline FHR (bpm)	110-150	100-110	<100			
		150-170	>170			
Variability						
Amplitude (bpm)	5-25	5-10 for >40 min >25	<5 for > 40 min			
Decelerations	None	Variable	Severe variable Severe repeated early Prolonged			
			Late or Sinusoidal			

3. Classification and management based on guidelines by ACOG 2010



4. Tracing classification and management based on NICE guidelines

		NICE Classification	
FHR Feature	Normal	Non-reassuring	Abnormal
	Reassuring		
Baseline FHR	110-160	161-180	<100
(bpm)			>180
Variability	>5	<5 for 30-90 minutes	<5 for > 90 min
Decelerations	None or early	Variable decelerations < 60 in depth <u>and</u> lasting <60 sec Present for > 90 min With >50% of contractions	Non-reassuring variable decels persist 30 min after starting conservative measures With > 50% of contractions
		Or Variable decelerations >60 in depth <u>and</u> lasting >60 sec Present for up to 30 min With >50% of contractions Or Late decelerations Present up to 30 min With > 50% of contractions	Or Late decelerations Present over 30 min Do not improve with conservative measures With > 50% of contractions Or Bradycardia or single prolonged decel > 3min

NICE Clinical management								
FHR Classification	Normal Reassuring	Non-reassuring and suggests need for conservative measures	Abnormal and indicates need for conservative measures AND further testing	Abnormal and indicates need for urgent intervention				
Definition	ALL FHR features are Normal/ Reassuring	1 Non-reassuring feature and 2 normal/reassuring features	1 abnormal or 2 non reassuring features	Bradycardia Prolonged deceleration> 3 min with baseline < 100				
Interpretation	Normal	Associated with increased risk of acidosis. Acidosis unlikely if acceleration are present	Likely to be associated with acidosis	Likely to be associated with acidosis or rapid development of fetal acidosis				
Management		A Consider underlying causes If baseline > 160-Check for fever Start conservative measures Inform coordinating midwife and obstetrician	B Management from A Offer to take fetal blood sample Take action sooner than 30 min if late decelerations are accompanied by tachycardia or reduced baseline variability	C Start A Urgently seek obstetric help Prepare for urgent birth Expedite birth if persistence > 9 min Reassess and discuss with the women if heart rate if heart rate recovers < 9 minutes				

5. Classification and management by Parer et al 2006.

G=Green, B=Blue, Y=Yellow, O=Orange, R=Red. Each of the modifiers "mild" "moderate" and "severe" was defined in numeric terms.

			Tra	cing (Classific	ation					
Moderate Varia	bility										
	Decele	rations	Recu	rrent V	ariable	Red	current	Late		Prolong	ed
	None	Early	Mild	Mod	Severe	Mild	Mod	Severe	Mild	Mod	Severe
Baseline	140116	Larry	IVIIIG	11100	Jevere	IVIIIG	Mou	Jevere	IVIIIG	Wiou	30,000
Tachycardia	В	В	В	Υ	0	Υ	Υ	0	Y	Υ	0
Normal	G	G	G	B	Y	В	Ϋ́	Y	Y	Ϋ́	0
Mild Bradycardia	Y	Y	Y	Y	0	Y	Υ	0	Y	Υ	0
Mod Bradycardia	Y	Ϋ́	'	•	0		0	0			0
Severe Bradycardia	0	0			0		O	0			0
gerere Bradyearaid											
Minimal Variab	ility										
	Decele	rations	Recu	rrent V	ariable	Red	current	Late		Prolong	ed
	None	Early	Mild	Mod	Severe	Mild	Mod	Severe	Mild	Mod	Severe
Baseline		-									
Tachycardia	В	Υ	Υ	0	0	0	0	R	О	0	0
, Normal	В	0	Υ	0	0	О	0	R	О	0	R
Mild Bradycardia	0	0	R	R	R	R	R	R	R	R	R
Mod Bradycardia	0	0			R		R	R			R
Severe Bradycardia	R	R			R			R			R
0 haant \/aviahili	. .										
Absent Variabili	ty										
	Decele	rations	Recu	rrent V	'ariable	Re	current	Late		Prolong	ed
	None	Early	Mild	Mod	Severe	Mild	Mod	Severe	Mild	Mod	Severe
Baseline											
Tachycardia	R	R	R	R	R	R	R	R	R	R	R
Normal	0	R	R	R	R	R	R	R	R	R	R
Mild Bradycardia	R	R	R	R	R	R	R	R	R	R	R
Mod Bradycardia	R	R			R		R	R			R
Severe Bradycardia	R	R			R			R			R
Sinusoidal	R		•								
Marked variability	Υ										

Clinical Management								
Color	Risk of acidemia	Risk of evolution	Action					
Green	0	Very low	None					
Blue	0	low	Conservative techniques and begin preparation					
Yellow	0	moderate	Conservative techniques and increased surveillance					
Orange	Borderline Acceptably low	High	Conservative techniques and prepare for urgent delivery					
Red	Unacceptably high	Not a consideration	Deliver					

6. Management algorithm for Category II tracings by Clark et al 2012.

Nine branches lead to one of three possible managements.

