

Appendix of EFM Classification Methods


1. EFM scoring system based on Krebs et al 1979

FHR Feature	Krebs Score		
	2	1	0
Baseline FHR (bpm)	120-160	100-119	< 100
		161-180	>180
Variability			
Amplitude (bpm)	6-25	3-5 or > 25	<3
Frequency oscillations/min	>6	3-6	<3
Accelerations	>4	1-4	0
Decelerations	None Early	Moderate Variable	Late Severe Variable Atypical variable

2. Classification based on FIGO methods in 1987

FHR Feature	FIGO Classification		
	Normal	Suspicious	Pathological
Baseline FHR (bpm)	110-150	100-110 150-170	<100 >170
Variability			
Amplitude (bpm)	5-25	5-10 for >40 min >25	<5 for > 40 min
Decelerations	None	Variable	Severe variable Severe repeated early Prolonged Late or Sinusoidal

3. Classification and management based on guidelines by ACOG 2010



4. Tracing classification and management based on NICE guidelines

NICE Classification			
FHR Feature	Normal Reassuring	Non-reassuring	Abnormal
Baseline FHR (bpm)	110-160	161-180	<100 >180
Variability	>5	<5 for 30-90 minutes	<5 for > 90 min
Decelerations	None or early	Variable decelerations < 60 in depth <u>and</u> lasting <60 sec Present for > 90 min With >50% of contractions Or Variable decelerations >60 in depth <u>and</u> lasting >60 sec Present for up to 30 min With >50% of contractions Or Late decelerations Present up to 30 min With > 50% of contractions	Non-reassuring variable decels persist 30 min after starting conservative measures With > 50% of contractions Or Late decelerations Present over 30 min Do not improve with conservative measures With > 50% of contractions Or Bradycardia or single prolonged decel > 3min

NICE Clinical management				
FHR Classification	Normal Reassuring	Non-reassuring and suggests need for conservative measures	Abnormal and indicates need for conservative measures AND further testing	Abnormal and indicates need for urgent intervention
Definition	ALL FHR features are Normal/ Reassuring	1 Non-reassuring feature and 2 normal/reassuring features	1 abnormal or 2 non reassuring features	Bradycardia Prolonged deceleration > 3 min with baseline < 100
Interpretation	Normal	Associated with increased risk of acidosis. Acidosis unlikely if acceleration are present	Likely to be associated with acidosis	Likely to be associated with acidosis or rapid development of fetal acidosis
Management		A Consider underlying causes If baseline > 160-Check for fever Start conservative measures Inform coordinating midwife and obstetrician	B Management from A Offer to take fetal blood sample Take action sooner than 30 min if late decelerations are accompanied by tachycardia or reduced baseline variability	C Start A Urgently seek obstetric help Prepare for urgent birth Expedite birth if persistence > 9 min Reassess and discuss with the women if heart rate if heart rate recovers < 9 minutes

5. Classification and management by Parer et al 2006.

G=Green, B=Blue, Y=Yellow, O=Orange, R=Red. Each of the modifiers “mild” “moderate” and “severe” was defined in numeric terms.

Tracing Classification												
Moderate Variability												
	Decelerations		Recurrent Variable			Recurrent Late			Prolonged			
	None	Early	Mild	Mod	Severe	Mild	Mod	Severe	Mild	Mod	Severe	
Baseline												
Tachycardia	B	B	B	Y	O	Y	Y	O	Y	Y	O	
Normal	G	G	G	B	Y	B	Y	Y	Y	Y	O	
Mild Bradycardia	Y	Y	Y	Y	O	Y	Y	O	Y	Y	O	
Mod Bradycardia	Y	Y			O		O	O			O	
Severe Bradycardia	O	O			O			O			O	
Minimal Variability												
	Decelerations		Recurrent Variable			Recurrent Late			Prolonged			
	None	Early	Mild	Mod	Severe	Mild	Mod	Severe	Mild	Mod	Severe	
Baseline												
Tachycardia	B	Y	Y	O	O	O	O	R	O	O	O	
Normal	B	O	Y	O	O	O	O	R	O	O	R	
Mild Bradycardia	O	O	R	R	R	R	R	R	R	R	R	
Mod Bradycardia	O	O			R		R	R			R	
Severe Bradycardia	R	R			R			R			R	
Absent Variability												
	Decelerations		Recurrent Variable			Recurrent Late			Prolonged			
	None	Early	Mild	Mod	Severe	Mild	Mod	Severe	Mild	Mod	Severe	
Baseline												
Tachycardia	R	R	R	R	R	R	R	R	R	R	R	
Normal	O	R	R	R	R	R	R	R	R	R	R	
Mild Bradycardia	R	R	R	R	R	R	R	R	R	R	R	
Mod Bradycardia	R	R			R		R	R			R	
Severe Bradycardia	R	R			R			R			R	
Sinusoidal	R											
Marked variability	Y											

Clinical Management			
Color	Risk of acidemia	Risk of evolution	Action
Green	0	Very low	None
Blue	0	low	Conservative techniques and begin preparation
Yellow	0	moderate	Conservative techniques and increased surveillance
Orange	Borderline Acceptably low	High	Conservative techniques and prepare for urgent delivery
Red	Unacceptably high	Not a consideration	Deliver

6. Management algorithm for Category II tracings by Clark et al 2012.

Nine branches lead to one of three possible managements.

