


Advanced Perinatal Systems with Decision Support Technologies

PeriCALM®

Supporting Clinicians with Real-time Pattern Recognition and Analysis at the Bedside

If you're still counting out Montevideo units and reading paper strips, you need PeriGen®.



PeriCALM provides decision support for nurses with real-time, automated assessments of labor progress and fetal heart rate patterns.

No other system offers the tools or the value:

- Clinical decision support tools
- Interfaces to your EMR
- Site License—one-time cost to your hospital

Pattern Recognition

PeriCALM analyzes the data, clinicians manage the care.

PeriGen creates innovative, evidence-based software that empowers clinicians to easily deliver safer patient care using decision support technologies. We offer affordable systems that customers can expand as their needs grow. Over 16,000 clinicians in leading hospitals and healthcare systems across North America depend on PeriGen smart systems for safer patient care.

700 Alexander Park, Suite 100, Princeton, NJ 08540 • www.perigen.com • 877.700.4755

Copyright © 2012 by PeriGen, Inc. All rights reserved

MKT-MAT-0051, rev. 01.01, 2012-05 Rx Only



PeriGen®

Empowering Clinicians. Advancing Perinatal Systems.




Advanced Perinatal Systems with Decision Support Technologies

PeriCALM[®]

Providing Better Situational Awareness through Decision Support

When tracings are analyzed and presented in detailed and long-term views, clinicians can stay ahead of the curve.

With the support of advanced pattern recognition data, nurses can take timely, clinical action to provide safer patient care.



Curve

Patterns

PeriCALM provides:

- Real-time, automated assessments of labor progress and fetal heart rate patterns
- Automated calculation of Montevideo units
- Pattern recognition based on NICHD guidelines, with confidence levels to support each analysis
- Both a long-term view and a detailed view of the tracing that appear simultaneously on the screen
- Objective, evidence-based analyses and statistical data at the bedside

PeriCALM analyzes the data, clinicians manage the care.

PeriGen creates innovative, evidence-based software that empowers clinicians to easily deliver safer patient care using decision support technologies. We offer affordable systems that customers can expand as their needs grow. Over 16,000 clinicians in leading hospitals and healthcare systems across North America depend on PeriGen smart systems for safer patient care.

700 Alexander Park, Suite 100, Princeton, NJ 08540 • www.perigen.com • 877.700.4755

Copyright © 2012 by PeriGen, Inc. All rights reserved

MKT-MAT-0052, rev. 01.01, 2012-05 Rx Only



PeriGen[®]

Empowering Clinicians. Advancing Perinatal Systems.



Advanced Perinatal Systems with Decision Support Technologies

PeriCALM®

The Most Clinically Advanced System on the Market

PeriGen's systems support clinicians across the OB continuum.


PeriGen helps 16,000 clinicians throughout North America provide safer care to their patients:

- Decision support at the bedside helps nurses focus on the patients who need them most
- Systems are user friendly
- Clinicians can view tracings and analysis at the bedside, central monitor or home computer
- Customers can expand the system with nursing and physician documentation, physician orders and protocol bundles as department needs dictate


Components of the PeriCALM suite:



PeriCALM® Curve™
Analysis of Labor Progression



PeriCALM® Patterns™
Analysis of Electronic Fetal Monitoring



PeriCALM® Shoulder Screen™
Assessment of Risk for Shoulder Dystocia
with Brachial Plexus Injury



PeriCALM® Contractility Index
Analysis of Contraction Rates

PeriCALM analyzes the data, clinicians manage the care.

PeriGen creates innovative, evidence-based software that empowers clinicians to easily deliver safer patient care using decision support technologies. We offer affordable systems that customers can expand as their needs grow. Over 16,000 clinicians in leading hospitals and healthcare systems across North America depend on PeriGen smart systems for safer patient care.

700 Alexander Park, Suite 100, Princeton, NJ 08540 • www.perigen.com • 877.700.4755

Copyright © 2012 by PeriGen, Inc. All rights reserved

MKT-MAT-0053, rev. 01.01, 2012-05 Rx Only

PeriCALM[®] Plus[™]

Full Documentation with Decision Support for Nurses and Physicians

Staying ahead of the curve is easy—with an evidence-based PeriGen[®] system behind you.

PeriCALM Plus builds upon the tracing management, surveillance and archiving functions of PeriCALM with full nursing and physician documentation. In addition to real-time, automated pattern recognition, PeriCALM Plus generates timely prompts, alerts and reminders to assist clinicians in providing safe patient care.

Advanced Support for the Continuum of Obstetric Care:

- Documentation and support for Prenatal, Triage, Antepartum, Intrapartum, Post-Partum and Well Baby
- OB-specific nursing and physician documentation allows clinicians to record the delivery of care in a patient chart specifically designed to match clinical workflow
- Decision support at the bedside provides clinicians with objective data analysis and the ability to see trends over time
- Modular system easily expands to include physician orders and protocol bundles as department needs require
- Standard interfaces including ADT, Pharmacy, Lab, and Anesthesia

